

Work Order - Urethane Partial Hand



Technical O&P Solutions

837 Sargent Ave. Winnipeg, MB
R3E 0C1 Canada

Please contact MOT Inc. before sending any
custom fabrication projects.

1-877-395-0081

fabrication@myrdalorthopedics.com

Customer Info:

Company: Address:	PO#: Patient: Sex, Age: Weight, Height: Diagnosis:
Contact: Phone #: Fax #: Email:	

Info / Options (please circle choices, and fill in blanks)

Affected Side: L - R

Urethane Durometer: 80D - 90D

Color: Light Tan - Medium Tan - Dark Tan - Black - Custom

Desired Style: Mildebrandt (Functional Opposition) - Medium Definition (Finger-like appearance) - High Definition (must provide life cast)

Strap Style: Dacron w/ Velcro - Dacron Buckle

Cast: Unmodified - Modified

If Unmodified: Please indicate any problematic areas. A check socket or modified cast is preferred.

For Mildebrandt style: the optimal situation is to provide a check socket with a foam buildup to give a rough idea of the desired shape.

For Medium Definition style: we ask for a tracing of the sound side, pictures, and an AP / ML sketch of desired hand position.

For High Definition style: we ask for a life casting of a hand in the desired position. Great results are often achieved by casting a sibling / relative.

If a "check hand" or life cast is not provided, please use this space to sketch the desired shape (AP/ML views), describe intended use and provide any information you feel may be helpful.